




Hemo bioscience™
The science in Blood Banking

READY-TO-USE QUALITY CONTROL FOR GEL CARDS



QUALITY CONTROL FOR C.A.T.
COLUMN AGGLUTINATION TECHNOLOGY

Hemo bioscience MQC-CAT is the **ONLY READY-TO-USE** product currently available that is intended for quality control testing of routine blood bank reagents used in manual column agglutination technology (CAT) blood typing and antibody screening systems.

The kit includes positive and negative red cells and plasma, providing a cost-effective method for quality control testing of routine CAT testing reagents and cards.

Evaluated with Grifols DG Gel® 8 and Ortho Clinical Diagnostics ID-MTS™ Gel systems.

PRODUCT NAME	DESCRIPTION	PRODUCT CODE	PRODUCT SIZE
MQC-CAT	Blood bank QC for CAT methods	H9350	RBC-POS 1 x 4 mL RBC-NEG 1 x 4 mL PLASMA-POS 1 x 10 mL PLASMA-NEG 1 x 10 mL

FEATURES & BENEFITS

- Unique ready to use product
- 63 day expiry dating
- Vials with color coded caps – easy to use
- Comprehensive worksheet included
- Red cell antigen profile provided – perform quality control of Rh antisera
- Compatible with Grifols DG Gel® 8 and Ortho Clinical Diagnostics ID-MTS™ Gel systems

HEMO BIOSCIENCE COMMITMENT TO SERVICE

CUSTOMER SERVICE

Hemo bioscience offers you flexible ordering options by fax, phone or email. To assist you in inventory management, we can set up standing orders or blanket orders. Most products are available for immediate overnight shipping. We develop, manufacture and inventory products in our facility, located in the heart of Research Triangle Park, North Carolina.

TECHNICAL SUPPORT

Hemo bioscience technical support is staffed by professional, qualified technologists who have donor center, immunohematology reference laboratory and transfusion laboratory experience; we understand your challenges in the laboratory. We have two fully operational immunohematology laboratories on site with qualified staff who can resolve product questions and guide you in deciding which products fit your needs.

QUALITY POLICY

Hemo bioscience strives to provide our customers with products and customer service that meet or exceed their expectations. We have established a comprehensive Quality Management System which provides a framework for measuring and improving our performance based on FDA CFR 21 part 820 standards, and we are ISO 13485 accredited. We have systems and procedures in place to support us in our goal of ensuring our customer requirements are met.



4022 Stirrup Creek Dr.
Suite 311
Durham, NC 27703

Tel: 1 866 332 2835 info@hemobioscience.com
Tel: 1 919 313 2888 www.hemobioscience.com
Fax: 1 919 313 2894

