

RARE ANTISERA (POLYCLONAL)

Hemo bioscience is the exclusive distributor of several rare polyclonal antisera. These FDA approved rare antisera are manufactured in the United States by the American Red Cross Diagnostic Manufacturing Division.

PRODUCT NAME	DESCRIPTION	CATALOG NUMBER	PRODUCT SIZE
Anti-Kell	IAGT	H500	1 x 5 ml
Anti-S	IAGT	H501	1 x 5 ml
Anti-Fy ^a	IAGT	H503	1 x 5 ml
Anti-M	Direct agglutinin	H504	1 x 5 ml
Anti-s (small)	IAGT	H505	1 x 5 ml
Anti-k (Cellano)	IAGT	H502	1 x 3 ml
Anti-Kp ^a	IAGT	H506	1 x 3 ml
Anti-Kp ^b	IAGT	H507	1 x 3 ml

FEATURES AND BENEFITS

- Standard indirect antiglobulin technique with 15 minute incubation
- Directly agglutinating anti-M
- Reliable supply of rare antisera from U.S. donors

HEMO BIOSCIENCE COMMITMENT TO SERVICE

CUSTOMER SERVICE

Hemo bioscience offers you flexible ordering options by fax, phone or email. To assist you in inventory management, we can set up standing orders or blanket orders. Most products are available for immediate overnight shipping. We develop, manufacture and inventory products in our facility, located in the heart of Research Triangle Park, North Carolina.

TECHNICAL SUPPORT

Hemo bioscience technical support is staffed by professional, qualified technologists who have donor center, immunohematology reference laboratory and transfusion laboratory experience; we understand your challenges in the laboratory. We have two fully operational immunohematology laboratories on site with qualified staff who can resolve product questions and guide you in deciding which products fit your needs.

QUALITY POLICY

Hemo bioscience strives to provide our customers with products and customer service that meet or exceed their expectations. We have established a comprehensive Quality Management System which provides a framework for measuring and improving our performance based on FDA CFR 21 part 820 standards, and we are ISO 13485 accredited. We have systems and procedures in place to support us in our goal of ensuring our customer requirements are met.



4022 Stirrup Creek Dr.
Suite 311
Durham, NC 27703

Tel: 1 866 332 2835
Tel: 1 919 313 2888
Fax: 1 919 313 2894

info@hemobioscience.com
www.hemobioscience.com

