



**Hemo bioscience™**  
*The science in Blood Banking*

## BLOOD BANK SPECIALIST SOLUTIONS

**Hemo bioscience offers a range of solutions for the treatment, testing, suspension and preservation of red cells.**

Our Alsever's can be used for in-house manufacture of screening cells, rare panel cells or any application that requires treatment, suspension and preservation of red cells. Alsever's red cell diluent contains no sodium azide, is salt adjusted for optimum red cell life, and helps preserve red cell antigen integrity.

Hemo bioscience's 0.2M DTT is a versatile general purpose reagent with numerous blood bank applications.

PRODUCT NAME	DESCRIPTION	CATALOG NUMBER	PRODUCT SIZE
Alsever's + ATP	Red cell diluent/ preservative	H339	1 x 1 L
0.2M DTT	Ready to use 0.2M Dithiothreitol	H411	1 x 2 mL

### FEATURES & BENEFITS

#### **Alsevers Solution**

- Traditional red cell diluent/preservative
- ATP added to extend life
- Convenient one liter pack

#### **0.2M Dithiothreitol (DTT)**

- Versatile general purpose reagent
- Applications include:
  - o Serum treatment to differentiate IgG and IgM antibodies
  - o RBC treatment to destroy selected blood group antigens

## HEMO BIOSCIENCE COMMITMENT TO SERVICE

### CUSTOMER SERVICE

Hemo bioscience offers you flexible ordering options by fax, phone or email. To assist you in inventory management, we can set up standing orders or blanket orders. Most products are available for immediate overnight shipping. We develop, manufacture and inventory products in our facility, located in the heart of Research Triangle Park, North Carolina.

### TECHNICAL SUPPORT

Hemo bioscience technical support is staffed by professional, qualified technologists who have donor center, immunohematology reference laboratory and transfusion laboratory experience; we understand your challenges in the laboratory. We have two fully operational immunohematology laboratories on site with qualified staff who can resolve product questions and guide you in deciding which products fit your needs.

### QUALITY POLICY

Hemo bioscience strives to provide our customers with products and customer service that meet or exceed their expectations. We have established a comprehensive Quality Management System which provides a framework for measuring and improving our performance based on FDA CFR 21 part 820 standards, and we are ISO 13485 accredited. We have systems and procedures in place to support us in our goal of ensuring our customer requirements are met.



**Hemo bioscience™**  
*The science in Blood Banking*

4022 Stirrup Creek Dr  
Suite 311  
Durham, NC 27703

Tel: 1 866 332 2835 info@hemobioscience.com  
Tel: 1 919 313 2888 www.hemobioscience.com  
Fax: 1 919 313 2894

