



**Hemo bioscience™**  
The science in Blood Banking

## VALIDATION AND COMPETENCY KIT



### VALIDATION AND COMPETENCY TESTING

PRODUCT NAME	DESCRIPTION	PRODUCT CODE	PRODUCT SIZE
Validation and Competency Kit	Samples for validation of automated and manual blood bank testing methods	H2010	10 x 5 mL

### FEATURES & BENEFITS

- Instrument ready with bar-coded label
- Suitable for tube, gel, and solid phase; manual or automated methods
- Includes positive and negative samples
- Range of clinically significant antibodies
- Variety of ABO types allows for validation of ABO reverse group testing
- Crossmatch compatibility testing validation using donor samples
- 5 mL fill volume – sufficient for multiple tests; method comparisons and competency testing
- 24 month expiry dating; open and closed vial

## HEMO BIOSCIENCE COMMITMENT TO SERVICE

### CUSTOMER SERVICE

**Hemo bioscience** offers you flexible ordering options by fax, phone or email. To assist you in inventory management, we can set up standing orders or blanket orders. Most products are available for immediate overnight shipping. We develop, manufacture and inventory products in our facility, located in the heart of Research Triangle Park, North Carolina.

### TECHNICAL SUPPORT

**Hemo bioscience** technical support is staffed by professional, qualified technologists who have donor center, immunohematology reference laboratory and transfusion laboratory experience; we understand your challenges in the laboratory. We have two fully operational immunohematology laboratories on site with qualified staff who can resolve product questions and guide you in deciding which products fit your needs.

### QUALITY POLICY

**Hemo bioscience** strives to provide our customers with products and customer service that meet or exceed their expectations. We have established a comprehensive Quality Management System which provides a framework for measuring and improving our performance based on FDA CFR 21 part 820 standards, and we are ISO 13485 accredited. We have systems and procedures in place to support us in our goal of ensuring our customer requirements are met.



4022 Stirrup Creek Dr.  
Suite 311  
Durham, NC 27703

Tel: 1 866 332 2835 info@hemobioscience.com  
Tel: 1 919 313 2888 www.hemobioscience.com  
Fax: 1 919 313 2894

